A stylized sunburst icon consisting of a central black circle with numerous thin, radiating lines extending outwards.

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Rx
ONLY

Foreword

This manual and the equipment it describes are for use only by qualified medical professionals who have been trained to use electrosurgical equipment.

Equipment covered in this manual:

Ultravision™ visual field clearing system, comprising:

| Item | Part number |
|---|-------------|
| Ultravision™ Generator | DAD-001-010 |
| Ultravision™ Battery (x2) | DSD-001-034 |
| Ultravision™ Battery Recharging Station | DAD-001-024 |
| Ultravision™ Mains Converter Power Supply* | DAD-004-009 |
| Ultravision™ Mains Converter* | DAD-004-012 |
| Ultravision™ Power Supply Unit | DSD-001-035 |
| Ultravision™ Patient Return Adaptor (Solid) 0.75M | DAD-001-006 |
| Ultravision™ Patient Return Adaptor (Split) 0.75M | DAD-001-007 |
| Ultravision™ Patient Return Adaptor (Solid) 1.5M | DAD-001-029 |
| Ultravision™ Patient Return Adaptor (Split) 1.5M | DAD-001-030 |
| Ultravision™ Patient Return Adaptor (Euro) | DAD-001-031 |
| Ultravision™ 5mm Trocar (single) | DAD-003-013 |
| Ultravision™ 5mm Trocar (x6) | DAD-003-014 |
| Ultravision™ 5mm Trocar (x60) | DAD-003-015 |
| Ultravision™ Ionwand™ SterilePack | DAD-001-003 |

*supplied separately or incorporated within a mains powered Ultravision™ visual field clearing system starter pack variant

Patents pending

Worldwide patents are pending for the Ultravision™ system, including PCT/GB2010/051196 and PCT/GB2012/052707.

Trademarks

Ultravision™ and Ionwand™ are trademarks of Alesi Surgical Limited.

Manufacturer

Alesi Surgical Limited Cardiff Medicentre Heath Park
Cardiff, CF14 4UJ, U.K.
<http://www.alesi-surgical.com>

For information call +44 (0) 29 2029 1022.

MEDICAL - GENERAL MEDICAL EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI
ES60601-1 (2005) + AMD 1 (2012), CAN/CSA-C22.2 No. 60601-1 (2014), IEC 60601-1-6 (2010, A1:2013) CAN/CSA-C22.2
No. 60601-1-6 (2011, A1:2015)

© Alesi Surgical Limited, 2020. All rights reserved.

Terms Used in this Manual

Warning

Indicates a situation which, if not avoided, could result in death or serious injury.

Caution

Indicates a situation which, if not avoided, may result in minor or moderate injury.

Be aware

Indicates that the operator should be aware of a hazard that may damage the product.

Advice

Indicates advice to the operator to optimise performance of the system.

Statement on Electromagnetic Compatibility

The Ultravision generator and battery recharging station require special precautions regarding Electromagnetic Compatibility (EMC) and need to be installed and put in to service according to the EMC information provided in Chapter 9.

Portable and mobile RF communications equipment can affect the Ultravision generator and battery recharging station.

Ultravision™ system warning

No modification of this equipment is allowed.

Warranty

Please consult term and conditions of sale for details of the warranty offered with the Ultravision system and its components.

Table of contents

| | |
|--|----|
| Foreword..... | 2 |
| Terms Used in this Manual | 3 |
| Statement on Electromagnetic Compatibility | 3 |
| Ultravision™ system warning | 3 |
| Warranty..... | 3 |
| Table of contents | 4 |
| Chapter 1: Ultravision™ – overview and general features | 6 |
| Indications for use. | 6 |
| Contraindication | 6 |
| Parts shipped with the Ultravision™ system | 6 |
| The Ultravision™ system – principle of operation | 7 |
| Operating Parameters | 8 |
| Transport and Storage | 8 |
| Patient, Operating Room and General Use Safety | 9 |
| Maintenance..... | 10 |
| Chapter 2: Controls, Sockets and Indicators | 11 |
| Ultravision™ generator – Front panel view | 11 |
| Ultravision™ generator – Side panel view | 12 |
| Ultravision™ generator – Rear view..... | 13 |
| Ultravision™ Battery Recharging Station Option – Topview..... | 13 |
| Ultravision™ Main Converter Option – Front panel view | 14 |
| Chapter 3: Preparing the Ultravision™ system for use..... | 15 |
| Periodic Inspection | 15 |
| Charging the battery using the Ultravision™ Battery Recharging Station..... | 15 |
| Inserting the battery into the Ultravision™ generator | 16 |
| Operation of the Ultravision™ Generator using the optional Ultravision™ Mains Converter..... | 16 |
| Inserting the Mains Converter into the Ultravision™ generator | 16 |
| Preparing Ultravision™ for use - Instruments requiring a patient returnelectrode | 17 |
| Preparing Ultravision™ for use – Instruments that do not require a patient return path | 17 |
| Chapter 4: Using the Ultravision™ system | 19 |
| Inserting the Catheter, trocar and Ionwand..... | 19 |
| Use of the Catheter, Trocar and Ionwand during surgery..... | 20 |
| End of surgery | 20 |
| Insertion of the Ultravision 5mm Trocar | 21 |
| Use of the Ultravision 5mm Trocar during surgery | 21 |
| End of surgical procedure..... | 21 |
| Chapter 5: Aftersurgery | 23 |
| Cleaning the Ultravision™ Patient ReturnAdaptor | 23 |
| Post surgery actions with the mains converter | 23 |
| Cleaning the Ultravision™ Generator | 23 |

| | |
|---|-----|
| Cleaning the Ultravision™ Battery Rechargingstation | 23 |
| Recharging the Ultravision™ battery..... | 24 |
| Disposal..... | 24 |
| Chapter 6:Troubleshooting..... | 25 |
| Chapter 7: Summary System Specifications | 27 |
| Chapter 8: Accessories and replacementitems | 29 |
| Chapter 9: Electromagnetic guidance | 30 |
| Symbols used on Ultravision | 178 |

Chapter 1: Ultravision™ – overview and general features

This section describes the indications for use for the Ultravision surgical smoke visual field clearing system and its features.

Caution

Read all warnings, cautions, and instructions provided with the Ultravision system before use. Failure to do so could result in death or serious injury to the patient or operator.

Read the following associated individual manuals as these contain safety relevant information.

| Part No. | Manual Description |
|-------------|--------------------------------|
| DLU-001-025 | Ionwand Pack IFU |
| DLU-001-026 | Patient Return Adaptor IFU |
| DLU-001-004 | Battery Recharging Station IFU |
| DLU-001-005 | Battery IFU |
| DLU-001-009 | Ultravision User Guide |
| DLU-003-004 | Ultravision 5mm Trocar IFU |
| DLU-004-001 | Mains Converter IFU |

Indications for use

The Ultravision system is indicated for the clearance of smoke and other particulate matter that is created during laparoscopic surgery.

Warning

Ultravision is not intended to be used in non-laparoscopic surgery procedures.

Do not allow the Ionwand to come in to contact with the uninsulated portion of an electro-surgical instrument during operation. In such a situation the Ultravision generator issues a continual audio-visual proximity indication. Continued contact for in excess of 10 minutes may result in an electrical discharge from the electro-surgical instrument that may be sufficient to cause atrial fibrillation depending on proximity to the thoracic cavity.

Use Ultravision with caution in the presence of internal or external pacemakers. Similar to electro-surgical devices, Ultravision can affect the performance of a pacemaker. Consult the pacemaker manufacturer or hospital Cardiology Department for further information when use of Ultravision is planned in patients with cardiac pacemakers.

Similar to electro-surgical devices, Ultravision may also interfere with the performance of an internal cardiac defibrillator (ICD). If the patient has an ICD contact the ICD manufacturer for instructions before performing an electro-surgical procedure with Ultravision.

Contraindication

None.

Parts shipped with the Ultravision™ system

The Ultravision system comprises both reusable and single-use, sterile disposable items. The single-use, sterile disposable items are supplied separately. Furthermore, the Mains Converter is an optional form of powering the generator that is supplied separately or within an alternate Ultravision system starter pack variant. Contact your local sales representative in order to establish a supply of these items. All other items required for Ultravision to function are shipped with the Ultravision system.

Before use, confirm that the shipping carton contains the following items:

- | | |
|---|----|
| 1. Generator | x1 |
| 2. Battery x 2 (battery powered system only) | x2 |
| 3. Recharging station (battery powered system only) | x1 |
| 4. Mains Power supply unit for recharging station (battery powered system only) | x1 |
| 5. Mains converter including power supply (mains powered system only) | x1 |
| 6. Solid-plate patient return adaptor | x1 |
| 7. Split-plate patient return adaptor | x1 |
| 8. User manual | x1 |

If any of the items listed above are not in the shipping carton, please contact your local sales representative before attempting to use the system.

Prior to use, inspect all of the above parts for any damage that may have occurred during shipping. If you have any concerns, please contact your local sales representative before using the system.

The following reusable items will require replacement:

| Item | Approximate number of uses before requiring replacement |
|--------------------------------|--|
| Battery | First to occur of 2 years, 200 charging cycles at 20°C (68°F) or 1 year from last recharge at 20°C (68°F). |
| Patient return adaptor (solid) | 500 mating cycles with Ultravision generator and electro-surgical generator |
| Patient return adaptor (split) | 500 mating cycles with Ultravision generator and electro-surgical generator |
| Battery recharging station | 1000 mating cycles with Ultravision battery |
| Mains Converter | 100 mating cycles with Ultravision generator |

Replacement reusable items are available from Alesi Surgical or its local distributors.

The Ultravision™ system – principle of operation

The Ultravision system removes the particulate material that is created by energy-based instruments during laparoscopic surgery from the visual field. It achieves this by electrostatic precipitation, a process that transiently charges the particulate matter and causes it to be rapidly precipitated inside the peritoneal cavity as it is created. The system has been specifically designed to improve the clarity of the visual field and to minimise the release of particulate matter and gases into the operating theatre environment during laparoscopic surgery.

The Ultravision system requires a patient return pathway, similar to monopolar diathermy but at 5 orders less current magnitude. The system has been designed such it can be configured in two ways, depending upon the type of cutting instrument that is used:


With monopolar instruments that require a patient return electrode

The Ultravision system is provided with a Patient Return Adaptor that allows the electro-surgical unit and the Ultravision generator to share a single patient return electrode with the monopolar instrument (Figure 1).

With bipolar and ultrasonic cutting instruments that do not require a patient return pathway

The Ultravision patient return socket in the generator allows a standard patient return electrode to be used and connected directly to the Ultravision generator (Figure 2).

The Ultravision system is not supplied with a patient return electrode but has been designed to be compatible with patient return electrodes that have an “international” style connector. Please contact Alesi Surgical or your local distributor if you experience any problems with your preferred patient return electrode.

The following components are required for Ultravision to function and are supplied with the system:

1. The generator unit

The generator unit is the source of the high voltage at low power. The output power is regulated automatically and cannot be adjusted manually. The unit is operated with power provided by a removable and rechargeable battery or an optional mains converter.

2 The battery

The generator is powered by a battery that allows approximately 10 hours of continuous generator operation. The system is supplied with two batteries and a recharging station. As with all batteries of this type, the life will depend upon the number of recharges. It is recommended that batteries should be replaced every 12 months or when a fully charged battery no longer delivers six hours of generator power.

3 The mains converter

Depending on customer preference, the generator can also be powered by a mains converter that allows the unit to be powered the mains supply using its bespoke power supply unit. This power supply unit has a choice of four different plug styles that allows it to be used in most countries and has a universal input allowing 90-264V AC. The Mains Converter is designed so that it can be stored within the battery bay of the Ultravision generator when not in use.

4 Battery recharging station

The Ultravision system is supplied with a twin-bay recharging station and a mains power supply unit different from that designed for use with the mains converter. The battery recharger mains power supply unit has four different plugs that allow it to be used in most countries and has a universal input allowing 90-264V AC.

5 Patient return adaptor – solid plate

The Ultravision system requires a patient return electrode. The patient return electrode can be shared with the solid or split plate return electrode of the electrosurgical generator if one is required, by using the patient return adaptor. The “solid plate” patient return adaptor is used with solid plate patient return electrodes. The blue end of the adaptor plugs into the socket on the front of the electrosurgical generator. The orange connector at the other end of the cable plugs into the front of the Ultravision generator. These are available with 0.75 or 1.5M leads.

6 Patient return adaptor – split plate

This is the same as the solid plate patient return adaptor but is used with split plates. It has the same central pin that activates the sensing function of the electrosurgical generator (if present). The split plate adaptor end is red in colour. These are available with 0.75 or 1.5M leads.

The following single-use, sterile disposable components are required for Ultravision to function and are supplied as separate items:

7 Ionwand™

The Ionwand is the active cable that carries the current from the Ultravision generator to the patient. It terminates in an atraumatic stainless steel brush that is responsible for the electrostatic charging of the particles created during surgery.

8 Catheter and trocar

The catheter allows the Ionwand to be inserted into the patient's peritoneal cavity during surgery. It is introduced percutaneously using the trocar.

9 Ultravision 5mm Trocar

The Ultravision 5mm Trocar is an alternative to the Ionwand™, catheter and trocar and comprises a laparoscopic port that can accept 5mm surgical instruments. This port also has a separate lumen that accepts a dedicated Ionwand™

Operating Parameters

The Ultravision system may be used within the following ranges

- Ambient temperature range 10°C to 40°C (50°F to 104°F)
- Relative humidity 20% to 75%, non-condensing operating room conditions
- Altitude below or equal to 3000m (above 70kPa)

Transport and Storage

- Ambient temperature range -10°C to 40°C (14°F to 104°F) except batteries
- Relative humidity 20% to 85%, non-condensing
- Battery ambient temperature range 0°C to 40°C (32°F to 104°F). Prolonged exposure to higher temperatures will significantly reduce battery shelf life.
- Battery shelf life – 1 year from last recharge at 20°C (68°F)
- Altitude below or equal to 3000m (above 70kPa)

Patient, Operating Room and General Use Safety

Warning

Only surgical staff trained to use electrosurgical equipment may operate the Ultravision system.

During initial set up and before each use, inspect the Ultravision system and battery recharging station for any damage. If damaged, do not use or attempt to repair. Call Alesi Surgical or your local distributor for assistance.

To maximise patient safety, contact between the lonwand and metallic instruments should be minimised. Otherwise, the metal instrument may be charged resulting in an electrostatic discharge to the patient or operator. This may be sufficient to cause atrial fibrillation depending upon proximity to the thoracic cavity.

If the Ultravision generator, Mains Converter, Power supply units or battery recharging station become wet for any reason, switch off and air-dry thoroughly before reusing.

Do not position the equipment so that it is difficult to disconnect the direct plug in Power Supply units from the mains.

Do not remove any covers or panels exposing the internal components of the Ultravision system and battery recharging station. Removal of covers or panels will void any Warranty.

Ensure that the battery is correctly inserted into the Ultravision generator. Failure to do so may result in a fire hazard. Ensure that the Ultravision™ Mains Converter is always inserted into the Ultravision Generator in the correct orientation. Failure to do so may damage the Mains Converter and equipment.

Do not use the Mains Converter to power equipment other than the Ultravision Generator. This may damage the Mains Converter and equipment.

Do not use any other third party Power Supply to power the Ultravision Generator. This may cause damage to the Ultravision generator and the third party power supply. The operator is also at risk of electrostatic shock or burn, due to charge build up in external capacitances.

Do not mechanically shock or drop a Mains Converter. This may damage or reduce the life of the Mains Converter.

Do not hold or suspend the Mains Converter or its power supply by its cable – this may damage the equipment.

Do not operate this device in potentially explosive environments, such as in the presence of flammable anaesthetics. The battery recharging station has only been certified for use outside the immediate operating theatre environment. It does not have sufficient redundant voltage or current protection to allow safe contact with the patient.

Use the battery recharging station and the mains converter with their specific mains power supply units. They are not mechanically interchangeable and provide different types of electrical protection.

Plug the power cable directly into the power socket without any extension cords or adaptor plugs.

The power cable assembly should be checked periodically for damaged insulation or connectors. Do not use a damaged power cable.

Do not connect a wet power cable to a power socket.

Inspect any sterile component for damage to the packaging before use. If any damage is present, do not use.

Do not reuse or re-sterilise the lonwand, trocar, catheter or the Ultravision 5mm Trocar. Doing so may result in cross contamination, infection or injury to the patient or medical staff and may damage the Ultravision system.

Caution

Read all warnings, cautions, and instructions provided with the Ultravision system before using. Do not stack equipment on top of the Ultravision system.

To maximise patient safety, contact between the lonwand and patient tissue should be minimised. Otherwise, performance

will be compromised and visibility will be impaired.

Be aware of trip hazards with the cable of the Ultravision Mains Converter.

Be Aware

For battery powered operation:

- Using any unapproved battery with the Ultravision system may cause damage and will void the warranty.
- Connect the battery charger power cable to a wall socket having the correct voltage. Otherwise, product damage to the recharger or battery may result.
- Power supply voltages below 90V AC will significantly reduce the rate of recharging.

For Mains Converter Operation:

- Power supply voltages below 90Vac (to the Mains Converter Power Supply unit) will cause unpredictable generator behavior (low battery warning or no operation).

Maintenance

Warning

Inspect the Ultravision generator, battery recharging station or mains converter before each use. If there is evidence of damage, do not use. Call Alesi Surgical or its local distributor for assistance.

Always turn off the Ultravision generator before cleaning.

Always unplug the battery recharging station before cleaning.

Always unplug the mains converter power supply unit and remove cable from the mains converter before cleaning.

Caution

Dispose of the Ionwand, trocar and catheter along with other operative waste materials, following the procedures for your institution.

Be Aware

Use only detergents that are already routinely used in your hospital to clean non-sterile electrosurgical generators and cables to clean the Ultravision generator, battery recharging station and patient return adaptors. Do not apply alcohol, caustic, corrosive, or abrasive cleaning or disinfectant compounds, solvents, or other materials that could cause damage.

Do not attempt to sterilise the Ultravision generator or patient return adaptors in any way.

Keep the Ultravision generator and battery recharging station away from liquid. These may damage internal components.


Chapter 2: Controls, Sockets and Indicators

This section describes all of the controls, sockets and indicators on the front and side panels of the Ultravision generator, Ultravision battery recharging station and Ultravision Mains Converter.

Warning

Read all instructions before operating this equipment. Failure to do so may result in damage to the generator and/or injury to the patient or operator.

Ultravision™ generator – Front panel view

- 1 – Power On/Off switch and indicator LED
- 2 – Patient return socket
- 3 – Ionwand socket
- 4 – Proximity indicator LEDs
- 5 – Battery level indicator
- 6 – Generator fault indicator LED

Warning

All third-party equipment used with Ultravision must be certified according to the relevant particular standards associated with IEC 60601. The user is responsible for verifying that any accessory equipment complies with the requirements of IEC 60601.

Be Aware

Connect accessories to the proper socket type. Otherwise, the system will not function properly.

1 – Power On/Off Switch and indicator LED

The power switch that controls power to the generator is located in the lower-left corner of the front panel (see Front Panel View, above). To switch the generator on, place the power switch in the ON position (I). The indicator light should illuminate green.

Turn the generator power off by placing the power switch in the OFF position (O). The indicator light should no longer be lit.

2 – Patient return socket and proximity indicator LED

The solid or split plate patient return electrode or Ultravision Patient Return Adaptor is inserted into this socket in order to establish the return path.

Be Aware

The Ultravision Patient Return Adaptor must be inserted the correct way up. Failure to do so may cause damage to the connector or Socket.

3 – Ionwand socket

The Ionwand connector is inserted into the Ionwand socket. This delivers the power to the Ionwand and causes the electrostatic precipitation of the particulates.

Be Aware

The Ultravision lonwand must be inserted in the correct orientation. Failure to do so may cause damage to the connector or socket.

4 – Proximity indicator LEDs

When the generator is switched on, the proximity indicator LEDs illuminate green. In order to function correctly, contact between the lonwand and patient tissue or other surgical instruments should be minimised. Both proximity indicator LEDs will flash when the lonwand comes within approximately 2mm (1/12”) of patient tissue. There is also an audible indicator that will sound. When the contact is removed the LEDs will stop flashing and the audible alarm will stop sounding.

Warning

Do not allow the lonwand to come in to contact with the uninsulated portion of an electrosurgical instrument during operation. Continued contact for in excess of 10 minutes may result in an electrical discharge from the electrosurgical instrument that may be sufficient to cause atrial fibrillation depending on proximity to the thoracic cavity.

5 – Battery level indicator

With a new, fully-charged battery, the battery level indicator displays the approximate power remaining as follows:

| | |
|------------------|------------------------|
| Green LED | >8hr power remaining |
| Amber LED | 2-8hr power remaining |
| Red LED | <2hr power remaining |
| Red flashing LED | <30min power remaining |

The red flashing LED is also accompanied by an audible indicator.

Where the Ultravision Mains Converter is used to power the Ultravision system, the Green LED is permanently illuminated on this indicator when the unit is switched on. If the Yellow or Red LED is illuminated, then the Mains Converter should be replaced.

Caution

The battery should not be used for a surgical procedure when the red flashing LED is present and should be re-charged. Failure to do so may result in loss of power during the procedure, leading to impaired visibility.

Be Aware

The battery may lose energy during transportation, storage and delivery particularly with prolonged exposure to higher ambient temperatures. It is recommended that the battery be fully charged before first use of the Ultravision system.

It is recommended that the battery be replaced once it is no longer able to deliver 6 hours of power when fully charged. Contact Alesi Surgical or its local distributor to purchase replacements.


6 – Generator fault indicator LED

This LED is not illuminated during use. If there is a problem within the generator, the LED will illuminate red and flash. This is also accompanied by an audible indicator.

Warning

Do not attempt to use the Ultravision generator when the generator fault indicator is lit. Contact Alesi Surgical or its local distributor to arrange for repair or replacement of the generator unit.

Ultravision™ generator – Side panel view

1 – Battery/Mains Converter insertion mechanism

1 – Battery/Mains Converter insertion slot and release mechanism

The battery or Mains Converter is inserted into the slot on the right-hand side of the Ultravision generator. It is subsequently


ejected from the Ultravision generator by pulling the ejection trigger toward the front of the generator.

Caution

The battery or Mains Converter must be inserted in the correct orientation. If not, damage to the battery, Mains Converter and/or generator may result and create a fire hazard.

Risk of fire, explosion or burns. Do not heat batteries above 60°C.

Ultravision™ generator – Rear view

1 – Volume control

1 – Volume control for audible indicators

The Ultravision generator has four audible indicators that accompany the visual indicators:

| Indicator | Tone | Description | Frequency | Visual indicator |
|-------------------|-------------------|--------------------------|------------|------------------|
| Low battery | High (1047 Hz) | A single 1 second tone | 2 seconds | Low battery LED |
| Proximity | Medium (523.3 Hz) | A single 0.3 second tone | 6 seconds | Proximity LEDs |
| Generator fault | Low (261.6 Hz) | A single 1 second tone | 2 seconds | Fault (!) LED |
| Power ON reminder | Low (261.6 Hz) | A single 1 second tone | 15 minutes | None |


The generator sounds the three different tones in sequence when switched on. The volume of these indicators may be adjusted by twisting the volume control on the rear of the generator. The volume is increased by turning the volume control clockwise, and reduced by turning anti-clockwise.

Warning

Ensure that the audible indicator can be heard within the operating theatre environment, and that the output LED indicators are visible to the operating personnel. Failure to observe indications so may result in:

- Capacitive charging of electrosurgical instruments and the risk of atrial fibrillation.
- Loss of function caused by power failure or a fault within the generator.

Ultravision™ Battery Recharging Station Option – Top view

1 – Power and fault indicator
 2 – Battery docking bay 1 and bay 2

- 3 – Battery charging status indicators
- 4 – Power supply socket

1 – Power/fault indicator

The power indicator lights green when power is being correctly supplied to the recharging station. If there is a fault within the recharging station it may still charge the batteries but it may be possible to over-charge them. This is indicated by the light changing colour from green to red.

2 – Battery docking bay 1 and 2

The recharging station is able to recharge up to two batteries at the same time. Batteries should be inserted such that the battery terminals engage with the terminals on the recharging station.

3 – Battery charging status indicators

When power is being supplied to the recharging station and a battery is inserted in the correct orientation the status indicator will illuminate. The indicator will be orange until the battery is fully charged. When the battery is fully charged the light will extinguish.

4 – Power supply socket

The battery recharging station is supplied with a power supply unit. The power supply unit has different plugs that can be attached to the power supply unit. Select the correct plug for your local power supply.


Be aware

If the power supply unit supplied with the recharging station does not have a suitable plug that works with your local power supply, please contact Alesi Surgical or your local distributor.

Only use the power supply unit provided with the battery recharging station. The use of other power supply units may damage batteries and recharging station.

When the power/fault indicator is red there is a fault with the recharger that may damage the batteries and the recharger should not be used. Contact Alesi Surgical to purchase a replacement recharger.

Ultravision™ Main Converter Option – Front panel view

- 1. Input socket for DC Power supply jack.
- 2. Generator battery contact (positive)

Chapter 3: Preparing the Ultravision™ system for use

This section describes how to prepare the Ultravision system for use, specifically:

- Periodic inspection
- Charging a battery using the battery charging station
- Inserting and removing a battery
- Preparing the mains converter for use.
- Preparing the Ultravision system for use:
 - Use with instruments requiring a patient return electrode
 - Use with instruments that do not require a patient return electrode

Warning

Read all warnings, cautions and instructions provided with the system and its accessories before use. Failure to do so could result in death serious injury to the patient or operator.

Periodic Inspection

The Ultravision system should be visually inspected at least once every year regardless of use. This inspection should include checks for:

- Damage to the generator and recharging station.
- Damage to the battery recharging station power cable and plug.
- Damage to the mains converter and its power supply including the cable and jack.
- Damage to the patient return adaptors cable or connectors.
- Proper mating between the patient return adaptors and the Ultravision generator and electrosurgical generator.

Warning

Do not use the Ultravision system with any damaged component. Damaged cables must be replaced before use. If there is any damage to the Ultravision generator or battery recharging station, contact Alesi Surgical or its local distributor.

Charging the battery using the Ultravision™ Battery Recharging Station

A battery should be recharged when purchased and whenever the Ultravision generator battery level indicator is red or flashing red. To fully charge a battery from zero to full power takes approximately 12 hours. A fully charged, new battery should provide at least 10 hours use of the Ultravision system.

1. Inspect the battery recharging station power cable for any signs of damage.
2. If there is no evidence of damage, attach the power cable to the recharging station.
3. Attach the correct face plate to the power cable plug and insert into a grounded, hospital-grade wall socket and confirm that the power/fault indicator light is illuminated green.
4. Insert a battery into a docking bay, ensuring that the battery terminals engage fully with the recharging station terminals. Confirm that the battery charging indicator is illuminated orange.
5. Allow the battery to charge until the charging indicator is extinguished. Remove the battery from the docking bay.
6. Unplug the power cable from the wall socket.

Batteries should be stored in a dry, clean environment when not in use.

Warning

Connect the battery charger power cable to a properly grounded, hospital-grade power socket of the correct voltage. Plug the power cable directly into the power socket without any extension cables and/or adaptors. The use of extension cables may increase the risk of a fire hazard.

The battery recharging station has been certified for use outside the immediate operating theatre environment. It does not have sufficient reserve voltage or current protection to allow safe contact with the patient.

Caution

Batteries must be recycled and not incinerated or disposed of in general hospital waste disposal. Do not attempt to directly connect the negative and positive terminals of a battery.

The battery contains acid and no attempt should be made to disassemble a battery.

Be aware

If the power supply units supplied with the recharging station does not have a suitable face plate that works with your local power supply, please contact Alesi Surgical or your local distributor.

The recharging station includes temperature-compensated characteristics and has been designed for optimal charging and life preservation of the batteries supplied with your Ultravision system. The use of other recharging stations even if mechanically compatible is likely to dramatically shorten the operating life of your Ultravision system batteries. Do not use batteries other than those supplied by Alesi Surgical with the recharging station.

Only use the power supply unit provided with the battery recharging station. The use of other power supply units may damage batteries and recharging station.

Do not mechanically shock or drop a battery. This may damage or reduce the life of the battery.

Do not use the battery to power equipment other than the Ultravision generator. This may damage the battery and equipment.

When the power indicator is red there is a fault with the recharger that may damage the batteries and the recharger should not be used.

It is recommended that the battery is removed from the generator if the generator is not to be used for some time.

Advice

To maximise battery life, batteries should be returned to the recharging station at the end of use, such as at the end of each day, and should be stored in as cool an environment as is possible. Leaving a battery uncharged over a weekend or holiday break is not advised.

Do not discharge battery completely. This will greatly shorten battery life.

Inserting the battery into the Ultravision™ generator

To insert a battery into the generator:

1. Locate the battery slot on the right hand side of the Ultravision generator.
2. Ensure that the battery bay is empty by pulling the battery ejection trigger towards the front of the generator.
3. Ensure that the battery is oriented with the battery terminals closest to the generator, with the arrow facing upwards and pointing to the generator.
4. Insert the battery until it engages with the contact terminals in the generator.
5. Switch on power to the generator by pressing the on switch located on the front panel.
6. Switch off the device until set up is complete.

Caution

Ensure that battery is always inserted into the Ultravision generator and battery recharging station in the correct orientation. Failure to do so may damage the battery and equipment and may cause a fire hazard.

Operation of the Ultravision™ Generator using the optional Ultravision™ Mains Converter**Inserting the Mains Converter into the Ultravision™ generator**

To insert the mains converter into the generator:

1. Plug the power supply of the Mains Converter into the appropriate mains socket (this may require the use of one of the regional adapters supplied with the Mains Converter). If the power supply units supplied with the Mains Converter does not have a suitable face plate that works with your local power supply, please contact Alesi Surgical or your local distributor.
2. Insert the jack at the end of the power supply unit into the Mains Converter inlet socket.
3. Locate the battery slot on the right hand side of the Ultravision generator.
4. Ensure that the battery bay is empty by pulling the battery ejection trigger towards the front of the generator.
5. Ensure that the Mains Converter is oriented with the terminals closest to the generator, with the arrow facing upwards and pointing to the generator.
6. Insert the Mains Converter until it engages with the contact terminals in the generator.
7. Turn on the mains power switch.
8. Turn the Ultravision generator fascia switch on when required and confirm that the LED on the generator fascia

is lit green and that the start up tone cycle is completed correctly. If this does not occur then consult the troubleshooting section (Chapter 6).

Warnings and Cautions

Ensure that the Ultravision™ Mains Converter is always inserted into the Ultravision Generator in the correct orientation. Failure to do so may damage the Mains Converter and equipment.

Do not use the Mains Converter to power equipment other than the Ultravision Generator. This may damage the Mains Converter and equipment.

Do not use any other third party Power Supply including the Battery Charging Station Power Supply to power the Ultravision Mains Converter. This may cause damage to the Ultravision generator and the third party power supply. The operator is also at risk of electrostatic shock or burn, due to charge build up in external capacitances.

Do not mechanically shock or drop a Mains Converter. This may damage or reduce the life of the Mains Converter.

Do not hold or suspend the Mains Converter or its power supply by its cable – this may damage the equipment.

Preparing Ultravision™ for use - Instruments requiring a patient return electrode

Before use, confirm that the following components are available:

- Ultravision generator with the mains converter or an adequately charged battery
 - Patient return electrode with “international” style connector: solid or split type
 - Correct Ultravision Patient Return Adaptor: solid or split
 - Electrosurgical generator unit compatible with the patient return electrode
- Set up the system as follows:
1. Connect the patient return electrode to the patient using standard surgical practice. Either solid or split variants may be used.
 2. Connect the patient return electrode to the patient return adaptor. Select the solid patient return adaptor for use with a solid return electrode. Select a split plate adaptor for use with a split plate.
 3. Connect the patient return adaptor end (BLUE – solid plate; RED – split plate) to the electrosurgical generator.
 4. Switch on the electrosurgical generator and ensure that the patient return path is recognised by the electrosurgical generator.
 5. Connect the other end of the patient return adaptor (ORANGE) to the Ultravision generator. The connector must be inserted in the correct orientation.

If the electrosurgical generator does not recognise the return electrode when using the patient return adaptor, follow the following procedure:

- Attach the patient return electrode directly into the electrosurgical generator. If this resolves the problem the patient return adaptor is worn and should be replaced.
- If this does not resolve the problem then replace the patient return electrode and repeat steps 1-4 above.
- If this does not resolve the problem then there is a problem with the electrosurgical generator.

Warnings and Cautions

Ensure that the selected patient return adaptor (solid or split) matches the patient return electrode that is attached to the patient. The patient return adaptor allows a solid patient return electrode to be inserted into a split patient return adaptor. Incorrect use may bypass any return pad sensing functionality of the electrosurgical generator.

Ensure that the Ultravision Patient Return Adaptor is used in the correct orientation.

Advice

If the surgical procedure may require both instrument types i.e. those that require a return pathway (e.g. monopolar) and those that do not (e.g. bipolar or ultrasonic) it is recommended that the system is set up in this way. This avoids the need to delay the procedure whilst reconfiguring Ultravision so that it can share the same patient return electrode.

The Ultravision system has been verified to have no effect on the split plate patient return contact monitoring systems employed in the monopolar generators listed on page 27.

Preparing Ultravision™ for use – Instruments that do not require a patient return path

Before use, confirm that the following components are available:

- Ultravision generator with the mains converter or an adequately charged battery

- Patient return electrode with “international” style connector (solid or split type) Set up the system as follows:
 1. Connect the patient return electrode to the patient using standard surgical practice. Either solid or split variants may be used.
 2. Connect the other end of the patient return electrode directly to the Ultravision generator by inserting it into the return socket on the face of the Ultravision generator. The connector must be inserted in the correct orientation.

Advice

Patient return adaptors should be made readily available in case the surgical procedure demands the use of an instrument that requires a patient return electrode.

Chapter 4: Using the Ultravision™ system

This section describes the use of the Ultravision system to clear the visual field of surgical smoke. This step requires the following consumable items:

- Ultravision catheter, and trocar
- Ionwand

Or

- Ultravision 5mm Trocar

These items are supplied separately.

The use of these accessories is described separately below:

Inserting the Catheter, trocar and Ionwand

The catheter is inserted using the trocar. The catheter and trocar are supplied assembled and ready for use.

The peritoneal cavity must be insufflated using standard practice before attempting to insert the catheter. Once insufflation has been achieved:

1. Open the Ionwand sterile packaging.
2. Remove the catheter/trocar assembly and the Ionwand from the packaging. Remove the protective sheaths from the trocar/catheter and Ionwand cable and check the integrity of the components.
3. Introduce the catheter by positioning at 90° to the abdominal wall and piercing with the catheter/trocar assembly
 - a. Always insert after the abdomen has been insufflated.
 - b. Always insert using direct vision to avoid iatrogenic injury to abdominal structures.
4. Remove the trocar from the catheter by holding the catheter and rotating the trocar anti-clockwise and dispose of safely.
5. Carefully insert the Ionwand down the catheter until it protrudes from the tip and can be visualised internally within the patient with the laparoscopic camera.
6. Secure the Ionwand in place by holding the catheter and rotating the hub of the Ionwand clockwise until it tightens. Do not overtighten:
7. Adjust the height of the catheter so that the Ionwand is free from contact with patient tissue and surgical instruments.
8. Connect the Ionwand connector to the active socket on the face of the Ultravision generator.

Warning

Single use only – the components have not been designed for re-sterilisation and re-use of this device poses a severe risk of cross-infection and/or device malfunction.

If the sterile packaging is open or damaged, do not use and discard safely.

If there is evidence of damage to any of the components, do not use and discard safely.

Always insert the catheter or trocar under direct vision once the abdominal cavity has been insufflated. Failure to do so may result in damage to internal tissue, blood vessels or organs.

Avoid major arteries and veins when inserting the catheter or trocar.

Never attempt to re-insert a partially or fully withdrawn trocar in the abdomen.

Do not switch on the Ultravision generator before inserting the Ionwand into the catheter.

As with electrosurgical instruments, contact between the Ionwand and uninsulated metallic instruments should be minimised to prevent electrostatic discharge to the operator or patient that may also transiently interfere with an ECG signal.

Repeated contact of the Ionwand with a video laparoscope may lead to temporary interference with the visual display and/or damage to the laparoscope.

Caution

The sterile components have been design specifically for use with this system. Only use items supplied by Alesi Surgical.

Advice

Insert the catheter in a position where it will not impair the surgical procedure and where it is least likely to come in to contact with other surgical instruments. Smoke clearing is most efficient when placed between the laparoscope and the site of surgery.

Ensure that the tip of the catheter is at least 25mm (1") from patient tissue to reduce the risk of accidentally activating the proximity indicator.

Cover the lumen of the catheter with a finger prior to inserting lonwand to minimise loss of pneumoperitoneum. Ensure that the lonwand brush is straight before insertion into the catheter.

Use of the Catheter, Trocar and lonwand during surgery

Before starting to cut or coagulate tissue and creating smoke:

1. Ensure that the lonwand is free of surgical debris and not touching patient tissue or a surgical instrument.
2. Switch on the Ultravision generator. If the battery status light indicator is flashing red and the battery alarm is sounding, replace the battery with a fully charged unit.
3. Perform surgical procedure.
4. In the event that the proximity indicator sounds when the lonwand™ is not touching patient tissue or a surgical instrument, perform the following steps:
5. Switch off the Ultravision generator.
6. Remove the lonwand™ from the catheter by holding the catheter and twisting the lonwand anti-clockwise, ensuring the catheter is not displaced.
7. Clean the lonwand™ and brush filaments by enclosing it within a sterile swab moistened with sterile saline and gently wiping it, moving the swab from the insulated portion of the lonwand towards the metal brush. Repeat until all contamination has been removed and the brush is clean.
8. Re-introduce the lonwand into the catheter and secure by holding the catheter and twisting the lonwand in a clockwise direction.
9. Switch on the Ultravision generator and confirm that the indicator condition has resolved and that smoke clearing has been resumed.
10. If the alarm continues to sound, remove and replace the lonwand.

Warning

To maximise patient safety, contact between the lonwand and metallic instruments should be minimised. Otherwise, the metal instrument may be charged resulting in an electrostatic discharge to the patient or operator. This may be sufficient to cause atrial fibrillation depending upon proximity to the thoracic cavity.

Accidental activation of a monopolar electrode in contact with the lonwand should be avoided. Provided the lonwand cable is connected to the Ultravision generator there is no risk of accidental burn or shock to the patient or surgical staff, but there is no proven redundancy in this protective insulation.

As with electrosurgical instruments, contact between the lonwand and uninsulated metallic instruments should be minimised to prevent electrostatic discharge to the operator or patient that may also transiently interfere with an ECG signal.

Repeated contact of the lonwand with a video laparoscope may lead to temporary interference with the visual display and/or damage to the laparoscope.

If the Ultravision catheter becomes disengaged from its position during use replace with a new device.

The system should not be used unless the user has been trained in its use.

End of surgery

When the surgical procedure has been completed or there is no longer any requirement for the Ultravision to be active:

1. Turn off the Ultravision generator.
2. Disconnect the lonwand cable from the generator.
3. Carefully remove the lonwand from the catheter

4. Remove the catheter and apply dressing to the wound site.
5. Dispose of lonwand cable and catheter safely.

Warning

Ensure that the Ultravision generator is switched off before removing the lonwand. Failure to do this may lead to electrical charging of conductive objects and create a risk of static shock to the operator.

Do not reuse or re-sterilise the lonwand, trocar or catheter. Doing so may result in cross contamination, infection or injury to the patient or medical staff, or damage the Ultravision system.

Dispose of the lonwand, trocar and catheter along with other operative waste materials following the procedures for your institution.

Advice

Removing the catheter as the last step in the procedure allows a surgeon to close the other surgical port sites without rapid loss of peritoneal pressure.

Insertion of the Ultravision 5mm Trocar

1. Open the sterile pack and remove the contents using sterile field technique.
2. Remove the protective sheath from the lonwand and check the integrity.
3. Insert the main obturator into the trocar and ensure correct assembly. Ensure the lonwand obturator is in the trocar and luer lock is secure.
4. Make an abdominal incision that is slightly smaller than the diameter of the trocar to avoid trocar instability. Insert the trocar under direct vision.
5. Apply light and continuous, but controlled downward pressure to the obturator to introduce the trocar through the abdominal wall. Ensure that the entire word ULTRAVISION is visible to the laparoscope inside the peritoneum.
6. When the trocar is in the desired position, remove the main obturator and the lonwand obturator.
7. Fully insert the lonwand and secure in place by tightening the luer lock.
8. Connect the lonwand connector to the Ultravision generator.

Use of the Ultravision 5mm Trocar during surgery

1. Ensure that the lonwand is free of surgical debris and not touching patient tissue or a surgical instrument.
2. Switch on the Ultravision generator.
3. Perform the surgical procedure.
4. In the event that the proximity indicator sounds when the lonwand is not touching patient tissue or a surgical instrument, perform the following steps.
5. Switch off the Ultravision generator.
6. Remove the lonwand from the trocar.
7. Clean the lonwand cable and brush filaments by enclosing it within a sterile swab moistened with sterile saline and gently wiping it, moving the swab from the insulated portion of the lonwand cable towards the metal brush. Repeat until all contamination has been removed.
8. Re-introduce the lonwand into the trocar.
9. Switch on the Ultravision generator and confirm that the indicator condition has resolved.
10. If the alarm continues to sound, replace the lonwand.
11. Switch off Ultravision generator once surgical cutting is complete.

End of surgical procedure

1. Ensure that the Ultravision generator is off.
2. Remove the lonwand from the trocar.
3. Remove the trocar from the abdomen. If port site haemostasis is not present, appropriate technique should be used to achieve haemostasis.
4. Apply dressing to the wound site.
5. Dispose of all items safely.

Warnings

The Ultravision 5mm Trocar has been designed specifically for use with the Ultravision System. Only use items supplied by Alesi Surgical Limited.

The Ultravision 5mm Trocar should only be used by surgeons with adequate training and familiarity with laparoscopic techniques.

If the pack is open or damaged, do not use and discard safely.

If there is evidence of damage to any of the components, do not use and discard safely.

Single use only – the components have not been designed for re-sterilisation and re-use of this device poses a severe risk of cross-infection and/or device malfunction.

Care should be exercised when introducing or removing sharp-edged laparoscopic instruments to lower the risk of inadvertent damage to the seal system.

Always insert the lonwand into the patient via the Ultravision 5mm Trocar before switching on the Ultravision generator.

Insert the trocar under direct vision and avoid iatrogenic injury to abdominal structures, major arteries and veins.

As with electrosurgical instruments, contact between the lonwand and uninsulated metallic instruments should be minimised to prevent electrostatic discharge to the operator or patient that may also transiently interfere with an ECG signal.

Repeated contact of the lonwand with a video laparoscope may lead to temporary interference with the visual display and/or damage to the laparoscope.

Ensure that the Ultravision generator is switched off before removing the lonwand.

Dispose of the Ultravision 5mm Trocar following the procedures for your institution.

The seal of the Ultravision 5mm Trocar can accommodate instruments with diameter of 5mm only. Use of instruments with diameter less than 5mm can result in loss of pneumoperitoneum.

Chapter 5: After surgery

This chapter describes how to clean the reusable components of the Ultravision system, how to recharge the battery when required and recommended actions if the mains converter is used in place of the batteries.

The reusable components should be cleaned with a mild detergent that is approved and routinely used in your hospital for cleaning surgical generators and non-sterile medical cables and connectors.

Cleaning the Ultravision™ Patient Return Adaptor

The patient return adaptors may require cleaning from time to time.

The cables and surfaces of the connectors should be wiped with a mild detergent. Cables and connectors should be allowed to air dry thoroughly before being reused.

Post surgery actions with the mains converter.

1. Switch off the Ultravision generator.
2. It is recommended that the Mains Converter is stored inside the Ultravision Generator until its next intended use.
3. The power supply jack should be removed by locating the barrel of the male socket slightly recessed in the battery bay and pulling gently to remove.
4. Disconnect the power supply unit from the mains supply. The cable itself should not be subjected to mechanical force as this may damage the unit.

It is recommended that the power supply and its cable are stored off the floor where the device will not be damaged.

Cleaning the Ultravision™ Generator

The generator may require cleaning from time to time:

1. Ensure that the generator is switched off ("O").
2. Eject the battery and remove or where a mains converter is in use, remove the power supply jack.
3. Ensure that all cables are detached from the generator.
4. Wipe all external surfaces thoroughly with a soft cloth that has been dampened with a mild detergent.
5. Ensure that the generator is dry before reusing.

Cleaning the Ultravision™ Battery Recharging station

The recharging station may require cleaning from time to time:

1. Ensure that the recharging station is unplugged from the power supply.
2. Ensure that the battery charging bays are empty.
3. Wipe the external surfaces thoroughly with a soft cloth that has been dampened with a mild detergent.
4. Ensure that the recharging station is dry before reconnecting the power supply.

Warning

Always disconnect the patient return adaptor from the electrosurgical generator and Ultravision generator before cleaning. Always turn off and unplug the Ultravision generator before cleaning. Allow the patient return adaptor and Ultravision generator to air dry thoroughly before reusing.

Do not attempt to sterilise the patient return adaptor or Ultravision generator.

Keep the patient return adaptor and Ultravision generator away from liquids. Liquids that enter the generator may damage internal components. Do not spray liquid cleaner into the Sockets of the Ultravision generator. Only apply appropriate cleaning liquids when applied to a cloth.

Be Aware

Do not attempt to clean the patient return adaptor or Ultravision generator with alcohol or caustic, corrosive or abrasive cleaning or disinfectant compounds, solvents, or any other material that could scratch or cause damage.

Recharging the Ultravision™ battery

Follow the procedure described in Chapter 2 or 3 to recharge the battery.

Disposal

All equipment with the exception of the disposables must be recycled at a WEEE-approved treatment facility. Do not incinerate or dispose of in general hospital waste disposal.

Chapter 6: Troubleshooting

This Chapter considers situations that may arise and suggests possible causes and recommended actions.

| Situation | Possible cause(s) | Recommended action |
|---|--|--|
| Generator does not power up when switched on | No battery power remaining | Replace with fully charged battery |
| | Faulty generator or mains converter | Contact manufacturer or local distributor |
| | Incorrect power supply being used with the mains converter. | Replace power supply with correct power supply. |
| Fully charged battery does not support generator operation for expected duration. | Battery at end of life through use or damage. | Replace with new battery |
| Battery not charging | Power cable not supplying power to recharging station | Confirm battery recharger power indicator is green. |
| | Battery not inserted correctly | Confirm that battery is inserted correctly with the terminals in contact with the terminals on the battery recharging station. |
| | Battery exhausted | Purchase replacement battery. |
| Battery recharger indicators (x3) not illuminated | No connection to wall adaptor | Check connection |
| | Wall adaptor not connected to an acceptable mains voltage supply | Check mains supply acceptable |
| | Battery recharger or wall adaptor damaged | Contact manufacturer or local distributor |
| Battery recharging lights fail to illuminate orange when battery still discharged | Battery not properly seated in battery recharger cradle | Reposition battery until recharging light changes to orange |
| | Battery recharger damaged | Contact manufacturer or local distributor |
| Battery recharging light(s) continuously illuminated orange after 12 hours of recharging | Battery damaged and continuously dissipating charging energy | Replace with new battery |
| | Battery recharger or wall adaptor damaged | Replace with new recharger |
| Battery recharging light(s) continuously illuminated red | Interference from other EMC sources | Switch Battery Recharger mains supply off. Ensure leads are not in contact or wrapped around those of another system. Reapply mains supply to the Battery recharger. |
| Electrosurgical generator does not recognise patient return electrode when using patient return adaptor <i>OR</i> No smoke clearing and electrosurgical instrument not working | Incorrect patient return adaptor | Ensure chosen solid or split adaptor matches the patient return electrode being used in the procedure. |
| | Defective patient return electrode | Connect return electrode directly to electrosurgical generator – if the return electrode is recognised it is NOT defective. |
| | Defective patient return adaptor | Perform previous two actions. If these do not resolve the issue, replace the patient return adaptor. |

| | | |
|---|--|--|
| No smoke clearing but electro-surgical instrument working | Ultravision generator switched off or has no power | Check generator is switched on and charged battery correctly inserted. Check that the mains converter is correctly inserted and that its correct power supply is being used. Check that the jack at the end of the power supply cable is correctly inserted into the mains converter and that the power is turned on. |
| | Ultravision system incorrectly configured | Ensure system configured correctly for the electro-surgical instrument being used. |
| | Defective lonwand | Confirm patient return path is correctly configured. Deliberately touch patient tissue with the lonwand. If the proximity alarm DOES NOT sound, replace the lonwand and repeat. |
| | Defective patient return path | If above action does not resolve the problem replace the patient return adaptor and patient return electrode. |
| | Interference from other EMC sources | Switch Ultravision off. Ensure leads are not in contact or wrapped around those of another system. Restart Ultravision generator. |
| Proximity alarm sounding frequently | lonwand is too close to patient tissue | Locate lonwand with laparoscope and confirm it is not too close to patient tissue. If too close, reposition the lonwand. |
| | lonwand is too close to instrument or port | Locate lonwand with laparoscope and confirm it is not too close to port or trocar. If too close, reposition the lonwand. |
| | lonwand partially contaminated with fluid or smoke particles | Remove the lonwand and wipe with sterile surgical wipe. Refer to "Use during surgery" in Chapter 4: Using the Ultravision™ system |
| | lonwand is defective | If above actions do not resolve the situation, replace the lonwand. |
| Proximity alarm sounding constantly | Rate of smoke production is too high to allow a smoke-precipitating electric field to be maintained. | Reduce the rate of smoke production by cutting more slowly. |
| | lonwand is heavily contaminated with smoke particles | Remove the lonwand and wipe with sterile surgical wipe. Refer to "Use during surgery" in Chapter 4: Using the Ultravision™ system. |
| | Catheter is contaminated with smoke particles and is providing a return path from the lonwand | Insert the catheter further into the abdominal cavity. |
| | lonwand is defective | If above actions do not resolve the situation, replace the lonwand. |
| Audible/visual indicator behavior not described by any of the above scenarios. | Interference from other EMC sources | Switch Ultravision off. Ensure leads are not in contact or wrapped around those of another system and restart the Ultravision system. Reposition the lonwand and cable if this does not resolve the situation. |


Chapter 7: Summary System Specifications

| Item | Specification |
|---|--|
| Storage and transport requirements | |
| Generator | -10°C to 40°C at maximum 85% relative humidity. Altitude below or equal to 3000m (above 70kPa) |
| Battery recharging station | -10°C to 40°C at maximum 85% relative humidity. Altitude below or equal to 3000m (above 70kPa) |
| Battery | -10°C to 40°C at maximum 85% relative humidity. Air Transport (IATA-DGR) UN No: 2800 Classification: Class 8 Proper Shipping Name: BATTERIES, WET, NON-SPILLABLE electric storage. |
| Mains converter | -10°C to 40°C at maximum 85% relative humidity. Altitude below or equal to 3000m (above 70kPa) |
| Patient return adaptor | -10°C to 40°C at maximum 85% relative humidity. Altitude below or equal to 3000m (above 70kPa) |
| Sterile consumables | -10°C to 40°C at maximum 85% relative humidity. Altitude below or equal to 3000m (above 70kPa) |
| Operating requirements | |
| Generator | +10°C to 40°C at 20-75% relative humidity, altitude below or equal to 3000m (above 70kPa) |
| Battery recharging station | |
| Battery | |
| Mains converter | |
| Patient return adaptor | +10°C to 40°C at 20-75% relative humidity |
| Sterile consumables | +10°C to 40°C at 20-75% relative humidity |
| Compatible patient return electrodes | Valleylab Adult Non-REM Polyhesive (Catalogue No. E7506) Valleylab Adult REM Polyhesive II (Catalogue No. E7507) Erbe MONOPlate™ Single Pad (Catalogue Number 20193-076) Erbe NESSY (Catalogue No. 20193-083 (-084 in USA)) |
| Compatible electrosurgical generators | Valleylab Force Fx Erbe VIO 300D Karl Storz Autocon II 400 Olympus ESG 100 |
| Battery capacity and operating life | 12 Volt, 2.0 Amp hours Estimated shelf life of 12 months from purchase or 200 charging cycles at 20C |
| Recharger input voltage range | 90-264V AC |
| Mains Converter Power Supply input voltage range | 100-240V AC(±10%) |
| Mains Converter Power Supply output voltage range | 13.4V DC ±5% |
| Maximum battery recharger current consumption | 0.5A at 230V AC |

| | |
|---|--|
| Mains Converter Power Supply current consumption (standby) | 230V AC : $\leq 0.1W$ |
| Battery run time (hours) | Approximately 10 hours |
| Number of battery recharges before replacement | 200 at 20°C |
| Generator, Mains Converter and battery re-charging station EMC compliance standards | EN 60601-1-2: 2014 (See Chapter 9) |
| Transit packaging rating: System | ISTA Class 2A |
| Transit packaging rating: Manisn Converter | ISTA Class 3A |
| Ingress protection rating for generator | 1litre saline poured onto top surface of generator (IEC 60601-2-2: 2009 5 th Edition, clause 201.11.6.3). |
| Patient Return Adaptor rating | Comply with IEC60601-2-2:2009 ed. 5th clause 201.15.101.2 and 201.15.101.3. |
| Generator maximum audible alarm (dBA) | 60dBA |
| Insulation type/shock prevention | Class 2 |
| Weight/Dimensions: Generator | 2.43kg; 395mm x 350mm x 100mm |
| Weight/Dimensions: Battery recharger | 175g; 140mm x 100mm x 50mm |
| Weight dimensions: Mains converter | 460g (exc power supply unit):135mm x64mm x 23mm |

Generator output characteristics

The graph below illustrates the relationship between voltage and current at the end of the lonwand cable.

The operating current is influenced by the distance of the lonwand from the patient return path (patient tissue). As the lonwand approaches the patient tissue the impedance (resistance) falls. This causes the current to increase and the voltage to fall. The generator monitors the current and an active feedback mechanism collapses the current as it approaches 10µA, which is the maximum DC current that can safely be applied to a patient. As a result the voltage falls below a level that is sufficient to cause electrostatic precipitation. This is accompanied the proximity indicator, which informs the operator that the lonwand is touching or very close to the patient tissue and requires repositioning.

Chapter 8: Accessories and replacement items

The following replacement items are available from Alesi Surgical or its local distributor:

| Item | Number per pack | Product code |
|---|-----------------|--------------|
| Ultravision™ Battery Recharging Station (recharging station, power supply unit and power cable) | 1 | DAD-001-024 |
| Ultravision™ Power Supply Unit (power supply unit and power cable) | 1 | DSD-001-035 |
| Ultravision™ Battery | 2 | DSD-001-034 |
| Ultravision™ Mains Converter | 1 | DAD-004-012 |
| Ultravision™ Patient Return Adaptor (solid) | 1 | DAD-001-006 |
| Ultravision™ Patient Return Adaptor (split) | 1 | DAD-001-007 |
| Ultravision™ Ionwand™ Sterile Pack | 10 | DAD-001-003 |
| Ultravision™ 5mm Trocar | 6 | DAD-003-014 |

Chapter 9: Electromagnetic guidance

Electromagnetic emissions – Ultravision Generator

With reference to the 4th Edition of IEC 60601-1-2:2014 The Ultravision generator is intended for use in the Special Environment of high powered medical electrical equipment associated with powered surgical devices. The customer or the user of the Ultravision generator should assure that it is only used in such an environment.

| Emissions test Compliance Electromagnetic environment – guidance | Compliance | Electromagnetic environment - guidance |
|--|------------|---|
| RF emissions CISPR 11 | Group 1 | The Ultravision generator uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. RF emissions |
| RF emissions CISPR 11 | Class A | The Ultravision generator is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies | |

Electromagnetic immunity – Ultravision Generator

The Ultravision generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision generator should assure that it is used in such an environment.

| IMMUNITY test | IEC 60601-1-2 test level | Compliance level | Electromagnetic environment – guidance |
|---|---|--|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact ± 15 kV air | ± 8 kV contact ± 15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/ burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines 100kHz repetition frequency | ± 2 kV for power supply lines ± 1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ±1 kV Differential Mode ±2 kV Common Mode | ±1 kV Differential Mode ±2 kV Common Mode | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % U_T ; 250/300 cycle | 0 % U_T ; 250/300 cycle | Mains power quality should be that of a typical commercial or hospital environment. If the user of the electrosurgical Electrosurgical generator requires continued operation during mains power interruptions, it is recommended that the Ultravision generator be powered from the identified battery accessory. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | 30 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Electromagnetic immunity – Ultravision Generator

The Ultravision generator is intended for use in the electromagnetic environment of a professional healthcare facility. The customer or the user of the Ultravision generator should assure that it is used in such an environment.

| IMMUNITY test | IEC 60601-1-2 TEST LEVEL | Compliance level | Electromagnetic environment – guidance |
|---|--|--------------------------|---|
| Conduct- ed RF IEC 61000- 4-6 Radiated RF IEC 61000-4-3 Including immunity to RF wireless comms equipment | 3 V _{rms} 150 kHz to 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz 3 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz | 3V 6V 3V/m | Portable and mobile RF communications equipment should be used no closer to any part of the Ultravision generator, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P}$ 80MHz to 800MHz $d = 2.33\sqrt{P}$ 800MHz to 2.7GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:  |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. | | | |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. | | | |
| ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ultravision generator is used exceeds the applicable RF compliance level above, the Ultravision generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ultravision generator. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m. | | | |

Recommended separation distances between portable and mobile RF communications equipment and the Ultravision Generator

The Ultravision generator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ultravision generator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ultravision generator as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter | Separation distance according to frequency of transmitter | | |
|---|---|--------------------|--------------------|
| | m | | |
| | 150 kHz to 80 MHz | 80 MHz to 800 MHz | 800 MHz to 2.7 GHz |
| W | $d = 1.17\sqrt{P}$ | $d = 1.17\sqrt{P}$ | $d = 2.33\sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.37 | 0.37 | 0.74 |
| 1 | 1.17 | 1.17 | 2.33 |
| 10 | 3.68 | 3.68 | 7.37 |
| 100 | 11.66 | 11.66 | 23.33 |

Electromagnetic emissions – Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision battery recharging station should assure that it is used in such an environment.

| Emissions test Compliance Electromagnetic environment – guidance | Compliance | Electromagnetic environment - guidance |
|--|------------|--|
| RF emissions CISPR 11 | Group 1 | The Ultravision battery recharging station uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. RF emissions |
| RF emissions CISPR 11 | Class A | The Ultravision battery recharging station is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies | |

Electromagnetic immunity – Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision battery recharging station should assure that it is used in such an environment.

| IMMUNITY test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
|---|---|--|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact ± 15 kV air | ± 8 kV contact ± 15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/ burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines 100kHz repetition frequency | ± 2 kV for power supply lines ± 1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % U_r ; 250/300 cycle | 0 % U_r ; 250/300 cycle | Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ultravision battery recharging station requires continued operation during power mains interruptions, it is recommended that the Ultravision battery recharging station be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | 30 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

NOTE U_r is the a.c. mains voltage prior to application of the test level.

Electromagnetic immunity – Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision battery recharging station should assure that it is used in such an environment.

| IMMUNITY test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
|--|--|-----------------------|--|
| <p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p> <p>Including immunity to RF wireless comms equipment</p> | <p>3 V_{rms} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz</p> | <p>3V</p> <p>3V/m</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the Ultravision battery recharging station, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P} \text{ 80MHz to 800MHz}$ $d = 2.33\sqrt{P} \text{ 800MHz to 2.7GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| <p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ultravision battery recharging station is used exceeds the applicable RF compliance level above, the Ultravision battery recharging station should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ultravision battery recharging station.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p> | | | |

Recommended separation distances between portable and mobile RF communications equipment and the Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ultravision battery recharging station can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ultravision battery recharging station as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter | Separation distance according to frequency of transmitter | | |
|---|---|--------------------|--------------------|
| | m | | |
| W | 150 kHz to 80 MHz | 80 MHz to 800 MHz | 800 MHz to 2.7 GHz |
| | $d = 1.17\sqrt{P}$ | $d = 1.17\sqrt{P}$ | $d = 2.33\sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.37 | 0.37 | 0.74 |
| 1 | 1.17 | 1.17 | 2.33 |
| 10 | 3.68 | 3.68 | 7.37 |

| | | | |
|------------|-------|-------|-------|
| 100 | 11.66 | 11.66 | 23.33 |
|------------|-------|-------|-------|

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Advena Ltd., Tower Business Centre