

OLYMPUS[®]

Your Vision, Our Future

Reprocessing Documentation, Workflow and Asset Management System

UNIFIA

Connect. Communicate. Decide.



UNIFIA

Introducing UNIFIA: A Reprocessing Documentation, Workflow and Asset Management System

As the complexity of cleaning and disinfecting endoscopes increases, staff compliance to endoscope reprocessing instructions becomes even more challenging. This requires staff supervisors, infection prevention champions and supply chain managers to have greater visibility into the utilization, cleaning, reprocessing and storage of the flexible endoscopes in their departments.

Connect


With UNIFIA, connecting and collecting data from reprocessing and storage devices has never been easier.

- UNIFIA's wireless scanner technology puts the ability to capture data from virtually all automated endoscope reprocessors, storage cabinets, leak testers and bio burden test kits in the palm of your hand.
- OER-Pro users can gain greater insight into their reprocessing times, scope serial numbers and HLD temperatures.

Communicate

UNIFIA communicates information from five locations, which include the reception desk, procedure room, storage locations, and soiled and reprocessing rooms.

- The UNIFIA Daily Dashboard and records are used to communicate critical patient and endoscope interactions across your entire GI department, keeping your staff in sync and working efficiently.
- A consolidated near-real-time Daily Dashboard view provides insight into critical operational, reprocessing and storage steps, including procedure room turnover, delayed reprocessing requirements, manual cleaning times and endoscope hang times.
- UNIFIA's wireless scanner technology sends onscreen alerts when scanning steps are missed, driving repeatable and compliant reprocessing processes while making sure all required data is collected.



Decide

The UNIFIA Analysis reporting engine provides users with the power to see their GI department like never before.

- Report templates offer vivid and easy to understand graphs and charts.
- Bring great transparency into your endoscope utilization and mix.
- Allow supply chain manager to make informed decision during endoscope acquisitions and service contract negotiations.
- Audit staff workloads and endoscope reprocessing compliance.
- Quick and easy endoscope reprocessing breach analysis.

Olympus is a registered trademark of Olympus Corporation, Olympus America Inc., and/or their affiliates.