

## Microsoft® Update and PHP Update Monthly Notification

**Summary:** This notification contains information about Microsoft security updates released this month and approved for implementation on the Image Stream Medical EasySuite®, LiveStream®, nCare®, and VaultStream® platforms. The only updates listed are those considered to be of a security or critical nature. Only the updates listed below have been validated with the EasySuite and EasySuite 4K, nCare MX/RX, EasyLink Router, LiveStream StatusBoard Engine, and VaultStream systems.

Also included is information on the latest PHP updates that are approved for the v8.x and up VaultStream systems.

**Products:** EasySuites and nCares, connected to a network and VaultStream, as well as the LiveStream StatusBoard Engine.

**Date:** August 12, 2020

### **Critical Security Bulletins (Obtained from Microsoft):**

**NOTE:** Image Stream Medical **does not approve** MS updates **KB975560, KB975562, KB2631813 or KB2780091** and strongly suggests they be excluded from automatic update scenarios. Failure to exclude these updates will result in equipment performance issues.

### **Approved Updates:**

August 2020 Updates		
KB Number	Platform	Title
<b>VaultStream v10.x Updates – Windows® Server 2012</b>		
KB4571703	VaultStream 2012	2020-08 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems
KB4570508	VaultStream 2012	2020-08 Security and Quality Rollup for .NET Framework 3.5, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 8.1 and Server 2012 R2 for x64
<b>VaultStream v10.x Updates – Windows® Server 2016</b>		
KB4569746	VaultStream 2016	2020-08 Cumulative Update for .NET Framework 4.8 for Windows Server 2016 for x64
KB4571694	VaultStream 2016	2020-08 Cumulative Update for Windows Server 2016 for x64-based Systems
KB4494175	VaultStream 2016	2020-01 Update for Windows Server 2016 for x64-based Systems
KB4565912	VaultStream 2016	2020-07 Servicing Stack Update for Windows Server 2016 for x64-based Systems
KB4565511	VaultStream 2016	2020-07 Cumulative Update for Windows Server 2016 for x64-based Systems
KB890830	VaultStream 2016	Windows Malicious Software Removal Tool x64 - v5.82
<b>nCare MX v10.2.2</b>		
KB4571694	nCare MX	2020-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4569746	nCare MX	2020-08 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
<b>EasySuite 4K v10.3   EasySuite 4K Controller (a.k.a. BlizCon)</b>		
KB4571694	EasySuite 4K Controller	2020-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4569746	EasySuite 4K Controller	2020-08 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
<b>EasySuite 4K v10.3   EasySuite 4K Video Manager (a.k.a. nCare RX2)</b>		
KB4571694	EasySuite 4K Controller	2020-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems

KB4569746	EasySuite 4K Controller	2020-08 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
<b>EasySuite 4K v10.3 – LiveStream StatusBoard 2</b>		
KB4571694	LiveStream StatusBoard 2	2020-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4569746	LiveStream StatusBoard 2	2020-08 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
<b>MedPresence Cart v10.3</b>		
KB4571694	MedPresence Cart	2020-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4569746	MedPresence Cart	2020-08 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
<b>PHP Updates</b>		
5.6.40 and 7.1.33 (7.3.21 available upon request for v10.2.2 VaultStreams only)		

For further questions concerning this notice, please call Technical Support at 1-800-848-9024



**Image Stream Medical**

One Monarch Drive Littleton, MA 01460 U.S.A.  
[www.olympusamerica.com/olympus-systems-integration](http://www.olympusamerica.com/olympus-systems-integration)

Image Stream Medical, Inc., is a wholly owned subsidiary of Olympus Corporation of the Americas